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With the ARTZT neuro Vision Charts for children you can train the input and cooperation of your perceptual systems of vision, balance and body awareness and thus improve your concentration. Perceptual input provides the brain with a great deal of important information about you and your environment and enables targeted and effective action and movement. Loss of information or interference can lead to less than optimal results. In addition, overloaded perceptual systems, such as tense eye muscles, can lead to headaches, concentration problems and learning difficulties. Training with the vision charts can be thought of as fitness training for the perceptual systems. Children with deficits in these systems may initially experience mild dizziness, nausea or headaches. Therefore, start slowly (2 x 5 minutes per day), possibly in a sitting position, and increase training intensity week by week. If symptoms persist, consult your doctor or therapist or contact a trained neuroathletics coach.

Most charts are hung on the wall at eye level for training. Choose a distance where you can easily see the symbols. If possible, and if you can see clearly, train without glasses.

Near/far Charts

Purpose: Training accommodation (focusing on objects at different distances)

Exercise: Use scissors to separate the small letter board (near board) from the large letter board (far board). Hang the far board on the wall and hold the near board in your hand. Take turns reading a letter from the near and far boards and say it out loud. Only switch to the other board when you can see the letter clearly. You can also do the exercise with each eye separately. It is best to use an eye patch to do this.

Vergence Charts (10 small exercise cards to cut out)

Purpose: To train gaze stabilization in the near and far range (divergence and convergence), thus improving the absorption and retention of learning content.

Convergence exercise „Squinting“: Hold the card at eye level in front of your eyes at a distance of 20-30 cm and create a superimposition/fusion of the two images by squinting. A third image will then appear in the center as a combination of the two images. Try to keep this image in focus, even if you move the card around.

Divergence exercise „Dreaming“: Hold the card in the center at eye level in front of your eyes at a distance of 20-30 cm and imagine that you are looking through the card (relax your eyes). A third image will appear in the center as a combination of the two images. Try to keep this image clear, even if you move the card around.

Lying Eight

Purpose: To improve control of head movements and cooperation between the visual system and motor perception, or to improve eye-hand coordination and hand motor skills.

Exercise 1: Copy the board and place it on the table. Using a coloured pencil, draw the figure eight in a fluid way, trying not to paint over the outer lines. Change speed, direction and hand. If you use different colours, you will end up with a real „work of art“.

Exercise 2: You will need the Motion Guidance Laser for this exercise. It is attached to your head and points straight ahead. Use the dot of light from the laser to follow the line of 8 on the board on the wall.

Emotions Saccade Chart

Purpose: To train empathy, precise eye movements and pinpoint fixation to improve concentration.

Exercise: With your head still, move your eyes from face to face (from left to right and from top to bottom) and imitate the facial expressions shown.

Eye Balance Chart

Purpose: To develop eye movement and head movement control, balance training, interaction between eye fixation and head movement.

Exercise 1 Eye movement: Keeping your head absolutely still, follow the lines on the board with your eyes. Gradually increase the speed of your eye movements. Use the ladybirds as a guide.

Exercise 2 Purposeful eye jumps: Keeping your head absolutely still, jump from ladybird to ladybird with your eyes. Increase the speed of your eye jumps over time. Can you see the dots on the wings of the ladybirds?

Exercise 3 Balance: Now move your head or the tip of your nose along the lines, keeping your eyes on the face in the middle of the picture.

Arrow Chart

Purpose: Faster processing of information, concentration and attention, interaction of sight-speech-orientation.

Exercise 1: Say out loud the direction of the arrows in the order shown (from left to right and from top to bottom). If there is a double arrow, say both directions. The quicker you do this, the better your systems will work together.

Exercise 2: Say out loud the direction of the arrows and at the same time point with one arm in the direction of the arrow. If there is a double arrow, both arms point in one direction. Increase your speed as soon as you get it right.

Exercise 3 - Pros: Same as exercise 2, but now the colours come into play:

Red: counter movement, with double arrow only one direction of your choice | **Blue:** leg movement | **Grey:** eye movement

Coordination Charts (4 charts with different levels of difficulty)

Purpose: These boards focus on eye-hand-foot coordination and the ability to translate an abstract image into movement. They stimulate memory and complex thinking.

Board 1: The triangles represent your hands. Right triangle - You clap your right hand against the side of your right thigh. Left triangle - You clap your left hand against the side of your left thigh. Triangle on the line - You clap both hands sideways against both thighs.

Table 2: The squares represent your feet. You stand on the toes of your right foot (square on the right), your left foot (square on the left) or both feet (square on the line).

Boards 3 and 4: Of course, you can also combine feet and hands and assign special tasks to the orientation of the triangles.

Word Type Chart

Purpose: This chart can be used to learn the parts of speech. The interspersed movement exercises also train the memory.

Exercise 1: For each capital letter you name a noun with that initial letter, for the lower case letters an adjective or verb. Use the arrows to move (turn left, turn right, squat, stretch, jump).

Exercise 2: Same as exercise 1, but now the small blue letters are adjectives and the pink letters are verbs.

Concentration Chart bd

Purpose: This board improves concentration and the ability to fix letters.

Exercise: Look at each letter and mark all the b's in yellow or all the d's in green. How many letters have you found?

Maths Coordination Chart

Purpose: This board can be used to practice and consolidate the basic arithmetic operations of addition and subtraction. The interspersed movement tasks also train memory skills. Explanation of symbols: Arrow pointing right = plus | Arrow pointing left = minus | Arrow pointing up = jump | Arrow pointing down = squat.

Exercise: You start with a starting number, e.g. 5, and a summand, e.g. 1. You stand in front of the board, count out loud and do the exercises: Stretch jump, $5+1=6$, stretch jump, $6+2=8$, $8-1=7$, squat,...

Number Chart

Purpose: To practice and consolidate basic arithmetic operations.

Exercise: The numbers on the board represent arithmetic results. You decide which operation to practice, e.g. addition. Which arithmetic problem results in 8? $2+6=8$. You can also assign different arithmetic operations to the different colours.

Number Search Chart

Purpose: to train the eye muscles, to change and fixate the gaze, thus improving orientation and the interaction between seeing, speaking and orientation.

Exercise: Look for the numbers in sequence (1 - 20 or 20 - 1) and always say them out loud when you find them. You can also look for the numbers with the help of an adult. Ideally, you should also take a torch and shine it on the number you have found. If you are doing the exercise in pairs, your partner can check that you have found the number. You can also stand very close to the board and touch the numbers as you find them.